







Diesel2p imaging system



- **[FOV]** Ultra-large field-of-view (FOV) optimized for multiphoton mesoscopy. Encompasses multiple brain areas within a single FOV.
- **[DUAL]** Two completely independent resonant-galvo scan engines. It's like having two mesoscopes. Temporally multiplex to see two places simultaneously, or have a dedicated mesoscale optogenetic pathway. Or do imaging and optogenetics activation in both channels.
- **[OBJ]** Compatible with conventional multiphoton objectives from Nikon, Zeiss, etc. Adapters for all major systems are available. The scan engine is infinity corrected and provides a 32 mm diameter beam at the back aperture, far beyond the requirements of conventional two-photon objectives.
- **[INF]** Infinity corrected and fully compatible with temporal multiplexing, lightbeads, Bessel beams, holography, and more.
- **[ROT]** Rotating air objective with full 360 degree range. Flexible to accommodate a large range of experiments, with the convenience of air immersion. No water interface needed. Plenty of working distance to make imaging easy, and navigate around complex instrumentation.
- **[WRLD]** Already in use in over labs in at least six different countries around the world. Real world proven. A robust workhorse system that stays aligned over time, and can be future ready, with a range of customizations possible. Friendly community sharing tips, tools, and more.

Multiscale random access:
Complex scan pattern +
Large field-of-view scanning



System specifications and features

Specification	Value	Notes
Scan engine	Dual, independent	<i>For multiplexing, optogenetics</i>
Field-of-view (native)	5 mm x 5 mm	<i>25 mm² in area</i>
Field-of-view (AO)	7 mm diameter	<i>38 mm² in area with adaptive optics (AO)</i>
Objective rotation range	360 degrees	<i>Any angle, including inverted</i>
Numerical Aperture (NA)	0.54	<i>Optimized for deep tissue imaging</i>
Working distance	8.0 mm	<i>No water needed</i>
Correction collar	Standard	<i>Optimize for your work</i>
Wavelength range	475 nm - 1100 nm	<i>Optimized for 910 nm and 1050 nm</i>
Field flatness	< 25 µm	<i>Flat imaging field</i>
Resonant scan axis	1500 µm @ 16,000 lines/s	<i>Wide resonant scanning stripes</i>
Resolution	Submicron lateral, 6-8 µm axial	<i>Even resolution across the FOV</i>
Typical imaging power	< 100 mW per channel	<i>Usually 50 - 70 mW</i>
Detection	GaAsP PMT or SIPM	<i>High sensitivity, robust</i>
Footprint	75 cm x 120 cm	<i>Height is adjustable</i>
Power requirements	Typically < 4 kW	<i>Depends on configuration</i>
Full peer-reviewed report	Click here	<i>Open access publication</i>

Large FOV example

6 mm x 5 mm FOV: 5 stripes
at 9.8 µm/pixel (2048 x 512 pixels)
12.5 frames/s with dual beams

Multi FOV example

Four 1.5 mm x 0.6 mm areas: 4 stripes
at 1.2 µm/pixel (2048 x 512)
15.6 frames/s with dual beams